FLEXBAND PLUS



SPEED

products.

Stitch integration eliminates the need for manual stitching to facilitate tissue passage or fixation of the matrix when indicated.

OPTIMIZED SUTURE RETENSION


Integrated whip-stitches create a robust suture/FLEXBAND interface.

SUTURE MANAGEMENT

Each FLEXBAND™ Plus utilizes a blue/white striped suture on one end and a solid white suture on the opposing end to aid in procedural suture management.



FLEXBAND PLUS



PRODUCT Number	PRODUCT DESCRIPTION	SIZE
41051	FLEXBAND Plus	0.3x8cm
41052	FLEXBAND Plus	0.3x16cm
41059	FLEXBAND Plus	0.3x32cm
41053	FLEXBAND Plus	0.5x8cm
41054	FLEXBAND Plus	0.5x16cm
41057	FLEXBAND Plus	0.5x32cm
41055	FLEXBAND Plus	0.7x8cm
41056	FLEXBAND Plus	0.7x16cm
HS332	Artelon Suture	20" suture loop w/needle

FLEXBAND™ Plus brings greater OR efficiency into our kinematically and biologically advanced matrix. With integrated sutures, the matrix is ready to pass or secure right out of the package. Differentiated suture tail colors aid in stitch management. And the integrated whip-stitch pattern optimizes suture retension within the matrix.

WWW.ARTELON.COM

