

FLEXBAND® TWIST


FROM ARTELON



THE NEXT PARADIGM **TWIST**

Reinforce and biologically restore injured ligaments and tendons with our next-generation FLEXBAND anchor system. FLEXBAND® TWIST enhances anchor placement precision, provides increased pullout strength, and offers efficient interoperative adaptability to soft tissue repairs and reconstructions.





CLEAN DRILLING. A reinforced drill guide with positive stop and metal lined sleeve promotes clean drilling and soft-tissue protection.



EFFICIENT ADAPTABILITY

PROCEDURAL VERSATILITY.

The Modification Cube provides a versatile intra-op selection of additional anchors, allowing you to choose either 3.85mm or 5.0mm fixation based on procedural and patient need.

OPTIMIZED FIXATION

CONTROL & STRENGTH.

Precisely deliver and set tension with eyelet insertion, then twist-in the anchor body for final fixation. This evolved insertion allows for enhanced FLEXBAND engagement, while increasing pullout strength 57% (compared to the tap-in system)³.

MORE MATRIX.

TWIST Kits come in two FLEXBAND length offerings, 12cm and 30cm, to cover an enhanced range of indications.

BIOLOGICAL RESTORATION

NOT JUST A BRACE, A RESTORATION. Artelon's FLEXBAND was uniquely designed as a non-inflammatory matrix that is capable of quickly integrating with regenerating connective tissue¹. After providing time-zero strength, FLEXBAND technology quickly integrates to share in tensile loading; protecting the healing construct and driving the mechanical stimulation required for optimal tissue regeneration and remodeling¹. Ultimately, the FLEXBAND matrix is replaced by highly organized type-1 collagen over 4-6 years².

References: (1) Lijjensten et al. J Mater Sci Mater Med. 2002 Apr;13(4):351-9 (2) Data on File (3) Data on File

FLEXBAND TWIST.12

3.85x17mm PEEK Eyelet Anchor
5.0x17mm PEEK Eyelet Anchor
0.5x12cm FLEXBAND Plus Dynamic Matrix
Cannulated Drill Bit
1.4mm Guide Wire
Tissue Protector

P/N TW012

x 2
x 1
x 1
x 1
x 2
x 1

FLEXBAND TWIST DRIVER 3.85

3.85x17mm PEEK Eyelet Anchor
with TWIST Driver

P/N TA385

x 1

FLEXBAND TWIST DRIVER 5.0

5.0x17mm PEEK Eyelet Anchor
with TWIST Driver

P/N TA500

x 1

FLEXBAND TWIST.30

3.85x17mm PEEK Eyelet Anchor
5.0x17mm PEEK Eyelet Anchor
0.5x30cm FLEXBAND Plus Dynamic Matrix
Cannulated Drill Bit
1.4mm Guide Wire
Tissue Protector

P/N TW030

x 3
x 1
x 1
x 1
x 2
x 1

ORDERING

